



South Asia

Choose certainty.
Add value.

**Test Report No.DEL/F(C/M)/13/000895 PART-2
Dated 15/03/2013**

Issued to : **KENT RO SYSTEMS LTD.**
A-2, Sector 59, NOIDA – 201 309

Attention : Mr. Dinesh Garg

Customer Ref. No. : NA

Sample Particulars

Sample Name : RO Water Purifier (Kent Supreme)

Sample Description : NA

Batch No. / Lot No. & Manufacturing Date, Best before Date : S/N- KR130303481

Sample Package & Quantity : Tradepack & 1 No.

Sampled by : Customer

Sampling Protocol : NA

Sample Condition : Sample received in good condition

Date of Sample Receipt : Received on 06/03/2013 4:57 pm

Date of Analysis : 07/03/2013

Date of Completion : 15/03/2013

Note: The submitted sample is not Drawn by the Laboratory.

Deepak Arora
Key Account Manager

Authorised By

Pankaj Jaiminy
Asst. Vice President

Note: The test report is electronically generated. Hence original signature is not required. For any technical query, please contact Dr. Sumant Nayak at sumant.nayak@tuv-sud.in and for any complaint please contact Rakhi Verma at Rakhi.Verma@tuv-sud.in .

Note: (1) General Terms & Conditions as mentioned overleaf, (2) The results relate only to the items tested, (3) the test report shall not be reproduced except in full without the written approval of the laboratory Laboratory: (4) For details of the accredited scope please contact the laboratory or visit www.nabl-india.org

TUV SUD South Asia Pvt. Ltd.
C-153/1, Phase - I,
Okhla Industrial Area
New Delhi-110020, New Delhi, India.

Phone :0091 11 30889611
Fax :0091 11 30889595
E-Mail :sumant.nayak@tuv-sud.in
www.tuv-sud.in

Regd. Office:
TUV SUD South Asia Pvt. Ltd.
Off Saki Vihar Road,
Saki naka, Andheri (East),
Mumbai-400093. India

TUV®



South Asia

Test Report No.DEL/F(C/M)/13/000895 PART-2
Dated 15/03/2013

S. No.	TEST PARAMETER	TEST RESULTS	TEST METHOD
A	Other Claims		
1	Drained water observed, % v/v	Nil	Volumetric measurement
2	Checking of Dual storage tanks in water purifier (one for filtered water and other for rejected water).	Present	Visual observation
3	RO water produced, % of feed water	50.0	Volumetric measurement
4	Checking of Automatic Operation of filtration process	Automatic system is in place which maintains the maximum level of water in the filtered water tank wherein it shuts off if the filtered water tank is full and starts whenever the water level falls in the tank	Operational observation

Remarks: The submitted Kent Supreme Mineral RO Water Purifier complies with the customer's claim with respect to above test parameters only.

Authorized By

Chemical

- END OF TEST REPORT-